



AXCEL NARROW

P335

FWD Powerbase



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at Sales@Meritshealth.uk

Features

- 22.5" wide front wheel drive chair is easy to maneuver in tight spots, narrow hallways and door frames
- Tight 90° turns
- Navigating over uneven terrain is easy with FWD 14" drive tires, maintaining contact.
- Independent hydraulic, high-performance suspension
- Large 9" rear casters to navigate obstacles and rough terrain
- Bright LED lights, included
- Tool-free, rear battery access
- Power recline adjustments, 85° - 160°
- Powered by MK batteries, rated up to 19 miles per charge

Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder & mount

Options

- ASL alternative drive
- Therafin positioning products
- Power articulating footplate
- Push handles
- Proform headrest with EZ-Clamp bracket

Specifications

Product#	P335
Weight Capacity	300 lbs.
Drivetrain	FWD
Seat Width	14 - 20"
Seat Depth	14 - 20"
Backrest Height	18 - 22"
Seat-to-Floor Height	17.5 - 19.5"
Seat Elevation (Optional)	12"
Overall Length	39.6"
Overall Width	22.5"
Overall Height	41.5"
Turning Radius	26" w/o footplate
Ground Clearance	3"
Total Weight	350-385 lbs.
Max Speed	5.2 mph
Battery Size (2 required)	MK 60 Ah, M34 Gel
Range Per Charge	up to 19 mi.
Motor	2-Pole/222W
Controller	R-NET 90A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14"x 3" hybrid gel
Rear Wheels	9" hybrid gel
Color Options	<div><div></div> Black Noir</div> <div><div></div> Dove Gray</div> <div><div></div> Sapphire Blue</div> <div><div></div> Ruby Red</div> <div><div></div> Camo Green</div>



Features MK batteries



merits Enhancing life through mobility.

©Merits Health Products UK Ltd
3 Haugh Road, Burntisland, Fife KY3 0BZ
www.MeritsHealth.uk | +44 7443 876057 | Sales@MeritsHealth.uk



rev. 0225