



AXCEL P335

FWD Powerbase



Enhancing Life Through Mobility.

Contact your Territory Manager or email us at Sales@MeritsHealth.uk

Features

- Independent hydraulic, high-performance suspension
- Large 229 mm casters to navigate obstacles and rough terrain
- Bright LED lights, included
- Tool-free, rear battery access
- Power recline adjustment 85° - 160°
- Power tilt adjustment 0° - 50°
- Adjustable 152.4 mm shear reduction recline
- Powered by MK batteries, rated up to 30.6 km per charge


Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder

Options

- 304.8 mm (12") elevating seat
- ASL alternative drive
- Therafin positioning products
- Power articulating footplate
- Push handles
- Proform headrest with EZ-Clamp bracket



Specifications

Product#	P335
Weight Capacity	136.1 kg (300 lbs)
Drive train	FWD
Seat Width	356–559 mm (14–22")
Seat Depth	356–559 mm (14–22")
Backrest Height	457–559 mm (18–22")
Seat-to-Floor Height	445–495 mm (17.5–19.5")
Seat Elevation (Optional)	305 mm (12")
Overall Length	1006 mm (39.6")
Overall Width	622 mm (24.5")
Overall Height	1054 mm (41.5")
Turning Radius	660 mm (26" w/o footplate)
Ground Clearance	76 mm (3")
Total Weight	158.8–174.6 kg (350–385 lbs)
Max Speed	9.7 km/h (6 mph)
Battery Size (2 required)	60Ah, Group 34 GEL
Range Per Charge	up to 30.6 km (19 mi)
Motor	2-Pole/222W
Controller	R-NET 120A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	356 mm x 76 mm (14" x 3" hybrid gel)
Rear Wheels	229 mm (9" hybrid gel)

- Color Options
- Black Noir
 - Dove Gray
 - Sapphire Blue
 - Ruby Red
 - Camo Green

Features MK batteries



merits Enhancing life through mobility.

©Merits Health Products UK Ltd
 Unit A4 Electra Park, Electric Avenue, Birmingham B6 7EB.
www.MeritsHealth.uk | 0121 368 9760 | Sales@MeritsHealth.uk



REV. 0326

