

merits

ULTRA P325

POWER TILT / LIFT / AFP



***Enhancing Life Through Mobility.***

Contact your Territory Manager or email us at [Sales@MeritsHealth.uk](mailto:Sales@MeritsHealth.uk)

## Features

- Center of gravity tilt 0° - 50° posterior
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 406.4 mm – 558.8 mm (16" - 22"); depth adjustable in 25.4 mm (1") increments
- 609.6 mm (24") base width and seat-to-floor height as low as 419.1 mm (16.5")
- Rated up to 30.6 km (19 mi) per charge

## Options

- 304.8 mm (12") elevating seat
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- Oxygen tank holder
- USB charger

## Models



**POWER TILT**



**POWER TILT & AFP**



**POWER TILT LIFT & AFP**

## Specifications

Product#	P325
Weight Capacity	136.1 kg (300 lbs)
Drivetrain	MWD
Maximum Speed	8.05 km/h (5 mph)
Per-Charge Range	Up to 30.6 km (19 mi)
Battery Size (2 required)	55Ah, 22NF AGM 51Ah, 22NF GEL
Ground Clearance	64 mm (2.5")
Turning Radius	660 mm (26")
Overall Length	914 mm (36")
Overall Width	610 mm (24")
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 419–495 mm (16.5"-19.5") Captain: 508–584 mm (20-23")
Seat Width	406–559 mm (16"-22") adjusted in 1" increments
Seat Depth	406–559 mm (16"-22") adjusted in 1" increments
Backrest Height	457–711 mm (18–28")
Total Weight	143.1–174.3 kg (315–384 lbs)
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	356 mm × 76 mm foam-filled (14" × 3")
Casters	152 mm × 51 mm foam-filled (6" × 2")
Color Options	<input type="checkbox"/> White <input checked="" type="checkbox"/> Red <input checked="" type="checkbox"/> Blue <input checked="" type="checkbox"/> Charcoal Gray

Optional MK batteries



**merits** Enhancing life through mobility.

©Merits Health Products UK Ltd  
Unit A4 Electra Park, Electric Avenue, Birmingham B6 7EB.  
www.MeritsHealth.uk | 0121 368 9760 | Sales@MeritsHealth.uk



REV. 0326

merits