



VELOCITY P325

MWD Powerbase



Enhancing Life Through Mobility.

Contact your Territory Manager or email us at Sales@MeritsHealth.uk

Features

- Center of gravity tilt 0° - 50° posterior
- Bright LED lights, included
- Shear reduction recline 90° - 150°
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 410 mm – 560 mm; depth adjustable in 25.4 mm (1") increments
- 609.6 mm base width and seat-to-floor height as low as 420 mm
- Rated up to 29 km per charge


Accessories

Based on chair configuration)

- USB charger
- Cup holder
- Oxygen tank holder
- Cane & crutch holder
- Foldable rear basket

Options

- 304.8 mm (12") telescoping power seat lift
- Joystick knob choices
- Trunk support pad
- Directional swingaway
- Captain, Pan or AVID powered seat
- Multiple footplate options



Specifications

Product#	P325
Weight Capacity	136 kg (300 lbs)
Drivetrain	MWD
Seat Width	410 mm – 560 mm (16" – 22")
Seat Depth	410 mm – 560 mm (16" – 22")
Backrest Height	580 mm (23")
Seat-to-Floor Height	420 mm – 500 mm (16.5" – 19.5")
Seat Elevation (Optional)	300 mm (12")
Overall Length	910 mm (36")
Overall Width	610 mm (24")
Overall Height	1020 mm (40")
Turning Radius	660 mm (26")
Ground Clearance	64 mm (2.5")
Total Weight	145 kg – 183 kg (320 – 404 lbs)
Max Speed	9.66 km/h (6 mph)
Battery Size (2 required)	55Ah, 22NF Gel, 55Ah, 22NF AGM
Range Per Charge	Up to 29 km (18 mi)
Motor	2-Pole / 275W
Controller	P&G VR2
Charger	6 Amp
Gradient	6°
Front Wheels	150 mm × 50 mm – solid (6" × 2")
Drive Wheels	360 mm × 80 mm – foam-filled (14" × 3")
Rear Wheels	150 mm × 50 mm – solid (6" × 2")

Color Options	Green
	Metallic Black
	Pink
	Metallic Royal Blue
	Candy Apple Red
	Platinum Silver

Optional MK batteries



merits Enhancing life through mobility.

©Merits Health Products UK Ltd
 Unit A4 Electra Park, Electric Avenue, Birmingham B6 7EB.
 www.MeritsHealth.uk | 0121 368 9760 | Sales@MeritsHealth.uk



REV. 0326

